

## Hypochlorite Method (Worm Bleaching)

- o Make fresh hypochlorite solution

<u>30 ml</u>	<u>40ml</u>
- 2.4 ml 5N NaOH	- 3.2 ml 5N NaOH
- 24.6 ml dH <sub>2</sub> O	- 32.8 ml dH <sub>2</sub> O
- 3 ml Bleach	- 4 ml Bleach

- o Wash worms off plates with M9 buffer
- o I usually wash off 2 or 4 plates per 15 ml conical tube, depending on amount of worms on plates
- o Let worms gravitate to bottom of tube and remove volume off the top until 4 ml total volume left
- o Add 2 ml of hypochlorite solution to each tube and mix well
- o Incubate on roller for 20 minutes at 25°C (30 minutes at 20°C)
- o Centrifuge worms for one minute at 4K with a clinical centrifuge
- o Discard supernatant and wash worms with 4 ml of M9 buffer
- o Centrifuge again for one minute at 4K and wash again with M9 buffer
- o Discard supernatant and resuspend worms in 1 ml M9 buffer and put into a 1.5 ml microcentrifuge tube
- o Incubate worms and M9 overnight on roller at 20°C
- o The next day plate L1s on seeded plates

### M9 Buffer (1 L)

6 g sodium phosphate dibasic (22mM Na<sub>2</sub>HPO<sub>4</sub>) or 12 g 7-hydrate

3 g potassium phosphate monobasic (22mM KH<sub>2</sub>PO<sub>4</sub>)

0.25 g magnesium sulfate 7-hydrate (1mM MgSO<sub>4</sub>) or 0.122 anhydrous

5 g sodium chloride (85mM NaCl)

pH to 6.5 with HCl and autoclave