

Preparation of RBC Ghosts for Proteomic Analysis

RBC Isolation:

1. Draw blood:
 - fill syringe with 5cc Anti-CoAg
 - draw 35 cc blood
2. Divide equally into two 50 mL flacon tubes (~20 mL each).
3. Add PBS to ~50 mL and mix gently by inversion.
4. Spin @ 2500 RPM, 18°C for 15 min.
 - SUP = platelets + plasma proteins
 - Interface = WBCs
 - Pellet = RBCs
5. Remove supernatant by aspiration.
6. Resuspend in PBS, repeat spin, and remove supernatant.
7. Repeat 4X.

RBC Lysis:

8. Pool RBCs into one tube.
9. Add 50 mL lysis buffer-pefabloc and vortex.
 - (Need ~ 10 volumes for efficient lysis)
10. Divide into two 50 mL corex tubes and spin @ 8K, 4°C for 20' on SS34.
 - (10 min prerun to cool centrifuge to 4°C)
 - ? spin at 11,000 rpm for 15 minutes at 4°C?
11. Remove supernatant without disturbing loose pellet – **SAVE SUP**
 - (tilt tube, remove sup, remove white clump)
12. Resuspend in 10 mL lysis buffer-pefabloc and repeat spin.

13. Resuspend pellet in lysis buffer (without pefabloc), aliquot in microfuge tubes, and spin 11,000 rpm for 15 minutes.

14. Store pellet at -80°C until processing for mass spectrometry.

MATERIALS

acid-citrate-dextrose anticoagulant

syringe with 5cc Anti-CoAg

50 mL flacon tubes

PBS (10 mM NaCl, 155 mM KCl, 10 mM glucose, 1 mM MgCl₂, 2.5 mM KHPO₄, pH 7.4)

PBS + 2.5 mM Pefabloc

Aspirator

Lysis buffer (5 mM NaHPO₄, 1 mM EDTA, pH 7.4)

Lysis buffer + 2.5 mM Pefabloc

two 50 mL corex tubes

SS34 rotor, sorvall centrifuge

microfuge tubes

The murine RBC ghost preparation was further solubilized in 0.1% SDS and 6M urea and digested with trypsin, chymotrypsin, or AspN plus 1/500 trypsin. Following digestion, the samples were boiled for 10 minutes to inactivate the enzymes and inhibitors leupeptin 1M , chymostatin 20 M , PMSF 0.1mM were added and the digests were combined. The resulting mixture was acidified with 10% phosphoric acid and separated on an Integral 100Q separation system (Applied Biosystems, Foster City, CA, USA) using cation-exchange chromatography with a 100mm x 4.6mm column containing 5mm particle Polysulfoethyl-A resin with 300 Å pore (PolyLC, Columbia, MD, USA). The buffers used were: buffer A, 20mM KH₂PO₄, 25% acetonitrile, pH 3.0; buffer B 20 mM KH₂PO₄, 350 mM KCl, 25% acetonitrile, pH 3.0. A flow rate of 700ml/min was run with a gradient of 0-30% buffer B over thirty minutes then 20-100% over 20 minutes. Thirty fractions were collected, each 90 seconds long.

Preparation of Human Triton Insoluble Ghosts

RBC ghosts were prepared by hypotonic lysis and extensive buffer wash.

Five hundred microliters of 2% Triton X-100 in PBS were added to 500 microliters of RBC ghosts for a final concentration of 1% Triton X-100.

They were vortexed and incubated on ice for 30 minutes,

then centrifuged at 100,000 x g for one hour in a Beckman ultracentrifuge at 4oC.

The supernatant was removed and the pellet was resuspended in 1 ml of 1M NaCl, 10 mM NaHPO₄, 5 mM EDTA, pH 7.4.

The sample was vortexed and spun again at 100,000 x g for one hour.

The salt-extracted, Triton insoluble protein pellet was solubilized with 0.1% SDS and stored at -80oC until analysis.

It was then solubilized further with 6 M urea, digested with trypsin, chymotrypsin, or AspN plus 1/500 trypsin and subsequently combined for strong cation exchange chromatography as above collecting 15 fractions.